



THE ONLY DEVICE
in the WORLD
that FEELS your BODY

For non-surgical BODY SHAPING



● MICROVIBRAZIONE COMPRESSIVA®

endoSPHÈRES® THERAPY


ENDOSPHERES THERAPY® is a brand-new treatment harnessing the power of COMPRESSIVE MICROVIBRATION®, patented system engineered for non-invasive skin rejuvenation, body shaping and physical therapy. The method has been studied by influential medical research centers, that have certified the efficacy of the treatment.

The treatment is delivered using a cylinder made up of special 55 silicone spheres positioned in a honeycomb-like pattern, that have a specific density and diameter.

SENSOR TECHNOLOGY

Our innovative SENSOR technology increases the efficacy of the Endosphères Therapy treatment. How does this work?


The sensors positioned inside the delivery handpiece detect tissue resistance and they determine the correct pressure required to treat each area optimally. This means it is easier to carry out **INDIVIDUAL BESPOKE TREATMENTS** and ensure best results.



5 ACTIONS


The device uses COMPRESSIVE MICROVIBRATION to transmit low frequency vibrations creating a pulsed rhythmic action on the tissue. The wave resulting from the vibrations has 5 therapeutic benefits on the tissue. It provides a well-rounded treatment. Initially it restores the normal state of tissues and microcirculation and in later stages works on reshaping and contouring the body.

pain
RELIEF




Compressive Microvibration acts on mechanoreceptors to reduce pain and tissue inflammation. It improves tissue oxygenation, reducing painful swelling.

VASCULAR
action




The honeycomb arrangement of the spheres stimulates the vascular system, affecting the blood flow in the tissues. The "vascular workout" of the tissue improves microcirculation.

LYMPHATIC
drainage




Compressive Microvibration uses the rotational movement of the spheres and the direction of the cylinder on the tissue to activate the lymphatic system to work correctly. This helps to eliminate excess fluids that have accumulated as a result of poor circulation.

body
CONTOURING



The oscillations of the spheres create a push pull action which remodels the tissue. The shape and contour of the body is improved after only a few sessions, as well as reducing unsightly orange peel effect visible on skin.

MUSCLE
toning



Compressive Microvibration technology tones the muscles as they actively resist the compression. The vibration stimulates the muscles and can improve tone in the treated areas.




▶ **PATENTED METHOD**
100% MADE IN ITALY

▶ **AUTOMATIC BOOST**
WITH POWER UP TO 355Hz

▶ **LAST GENERATION**
TECHNOLOGY

▶ **PRESSURE**
SENSOR


▶ an **ADDITIONAL CYLINDER**
may be requested to treat specific
vascular conditions



35 shore silicone



15 shore gel



The low frequency mechanical movement of the spheres work deep into the tissue and act on the main causes of cellulite

▶ **ENSURES**
customer retention

AK SENSOR

▶ **PAINLESS**

▶ **100%**
non-invasive

▶ **VISIBLE RESULTS**
noticed after only
a few sessions

BOOST 355 Hz

DYNAMIC frequency/pressure
INDICATOR

REMOTE CONTROL JOYSTICK

BODY TREATMENT

Areas of application

Cellulite • Orange Peel skin • Muscle Contractures
Oedema • Swollen legs • Arms • Buttocks • Ankles • Belly fat
Breast tightening • Post-pregnancy treatment

- ▶ REDUCES LYMPHATIC FLUID RETENTION
- ▶ REDUCES ORANGE PEEL SKIN
- ▶ IMPROVES VENOUS CIRCULATION
- ▶ IMPROVES BODY SHAPE
- ▶ TONES MUSCLES

▶ **INCREASE OF**
salon/clinic revenue

▶ **EXCELLENT**
return on investment



▶ The FIRST device in the world to AUTOMATICALLY change the intensity of the microvibrations according to tissue resistance.

● MICROVIBRAZIONE COMPRESSIVA®

endoSPHÈRES® THERAPY

Clinical observations have confirmed **Endosphères Therapy** as a **valid and trustworthy treatment** providing a solution for cellulite, poor venous circulation, fluid retention and muscle tension.

In the studies, in order to obtain the correct diagnosis, patients were subjected to a clinical evaluation. This included the inspection and palpation of the tissue. Echo colour Doppler, echography, thermography and biometrical impedance analysis were carried out. Analysis was carried out after every six sessions to give an insight into the tissue changes. The results show the efficacy of **Compressive Microvibration**.



Courtesy
"G. d'Annunzio"
UNIVERSITY
CHIETI/PESCARA - ITALY
Dept. Medicine Physical
Therapy and Rehabilitation

CLINICALLY PROVEN RESULTS

Results after **12** sessions



Pictures: Archive of XI SIES INTERNATIONAL CONGRESS OF AESTHETICAL MEDICINE AND SURGERY

Ak
sensor



Ak Sensor

TECHNICAL DATA

Diameter	315mm(up)/364mm(down)
Height	1.150 mm
Weight	17 kg
Vibration	29 - 355 Hz
Cooling	ED/AIRCOOLING
Sensor System	FND-21K
Touch Screen	Capacitivo 7 Colour
Sound Level	56dBA
Electrical Features	da 100 a 240 V./ 50-60Hz

www.endospheres.com

The Endosphères Therapy Ak Sensor patent
belongs to the Fenix Group and is manufactured by Fenix Group



The Endosphères Therapy brand belongs to Fenix Group, a 100% Italian company. The group has been at the forefront in the field of aesthetics and medical-aesthetic equipment, always basing its strength on Research and development.

The Fenix Group has developed Endosphères Therapy, a new technology, using a patented Compressive Microvibration method. In recent years this has proved itself as the treatment of excellence to combat edematous fibrosclerotic panniculopathy (cellulite), adipose accumulations and face rejuvenation.

The Fenix Group is an innovative company investing in

marketing and PR campaigns both in trade press and consumer press. It is active across all social media platforms. Customers are assisted with in salon/clinic marketing to promote the treatments locally in their area. Service and technical support is available throughout the world.

The Fenix Group's mission is to design high quality electro-medical devices and manufacture the most innovative, effective and science-based technologies. These essential features ensure safety and efficacy making the Fenix Group a successful and customer-driven company.



Headquarters ITALY

Fenix Group s.r.l. - Strada Lungofino, 187 - c/o C.C. Ibisco B9/10 Città Sant'Angelo (PE)
Tel: +39 085 4483235 - info@fenixgroup.it



Branch Office U.K.

Fenix Group LTD - 51 Holland Street, London W8 7JB
Tel: +44(0)2070960040 - ukmanagement@fenixgroup.it



Branch Office SPAIN

Fenix Endospheres Therapy SL - Avda. Barón de Cárcer, 50/1 - Valencia
Tel: +34 627 508380 - esmanagement@fenixgroup.it



Branch Office FRANCE

Fenix France S.A.S. - Startway 43-47 Avenue de la Grande Armee - 75016 PARIS
Tel: +34(0)170395146 - frmanagement@fenixgroup.it

For other countries look for distributors on www.endospheres.com