


the latest anti-aging procedure





M.A.S.

sensorised ablative microvibration



For proper treatment of the face, the first step is exfoliation. The superficial horny tissue is thus purged of dead cells (the ones without nucleus).


Sensorized Ablative Microvibration, given by the inclination of the spheres, activates a friction mechanism on the epidermis that **smoothes** it and **facilitates its oxygenation**; it is most commonly known as **mechanical peeling**.

The skin is thus prepared for the next phase.

RESULTS*

Texture of the skin improved | **Tightening** of face oval Soft peeling | Increased **brightness**

EXFOLIATING EFFECT



mechanical peeling



15 shore gel

FACE


sensorized compressive microvibration

The FACE handpiece is the real innovation in the Eva by Endospheres method. Equipped with:


- a sensitive system of sensors that facilitate the operator's work by guiding him through LED lights positioned on the back of the handpiece. The correct pressure to be applied on each individual area is thus indicated according to the tissue detected;
- **gel microspheres** made with an exclusive, patented gel that allows it to act in depth while remaining gentle on the dermis for painless treatment:
- a self-lubricating sliding system, a mechanism devised by Fenix Group engineers, which allows the device to work without additional lubricants, uses the device without barriers between the skin and the spheres, amplifying the result.



ACTS ON MICRO-CIRCULATION



EFFECT ON
MUSCULAR
HYPERCONTRACTION



LYMPHATIC DRAINAGE



RESULTS*

Wrinkle reduction


Reduction in facial expression wrinkles

Tissue vascularization

Increased tissue and muscular tone

Improvement of skin and **muscle structure**






stimulated and reactive to receiving the nutrients that will be delivered. Facial muscles are turgid thanks to vascularisation. The epidermis is purified and pores

draws the substance into the tissue.

This technology facilitates the absorption of substances not only by superficial tissue, and reduces operators' working time.

CARRIER ACTION



*The results of the EVA by Endospheres are temporary and require maintenance treatments.

RESULTS*

Wrinkle reduction

Increase in thickness of tissue

Reduction in facial expression wrinkles



Los beneficios de la patentada Microvibración Compresiva se han adaptado al dispositivo Eva, **creación** de **Endospheres'** enteramente dedicada al cuidado facial. Here, the Aquí, la filosofía Endospheres se materializa en tres cabezales diferentes, uno para cada fase. La revolución en el tratamiento facial radica en la consecuencialidad de las funciones de cada cabezal, diseñadas para lograr resultados inmediatamente visibles. Los tejidos del rostro son mucho más sensibles que los de otras zonas del cuerpo, por lo que se han desarrollado microesferas de gel patentadas para tratarlos con suavidad y eficacia.

100% NO INVASIVO | REPETIBLE | EFECTIVIDAD DESDE LAS PRIMERAS SESIONES*



phase 1

M.A.S. sensorised ablative microvibration exfoliating effect



phase 2

FACE

sensorised compressove microvibration

acts on micro-circulation, lympatic drainage, muscle toning



phase 3

V.E.V.

vibro-electro-delivery

carrier action



AREAS OF APPLICATIONS

Lip area • Periocular area (crows feet) Forehead area • Genal-nasal area

Cheekbone area • Glabella area • Treatment on breast

Treatment on neck and neckline



Headquarters ITALY

Fenix Group s.r.l. - Strada Lungofino, 187 - c/o C.C. Ibisco B9/10 Città Sant'Angelo (PE) Tel. +39 085 4483235 - info@fenixgroup.it

Branch Office SPAIN

Fenix Endospheres Therapy SL - Valencia - espmanagement@fenixgroup.it

Branch Office U.K.

Fenix Group LTD - London - ukmanagement@fenixgroup.it

Branch Office FRANCE

Fenix France S.A.S. - Paris - frmanagement@fenixgroup.it

Branch Office SWITZERLAND

Fenix Group s.r.l. Switzerland - Bussigny - info.swiss@endospheres.com

Branch Office Deutschland

Fenix Group s.r.l. Deutschland- Stuttgart - germanagement@fenixgroup.it

Branch Office Benelux

Endospheres - Bruxelles - benelux@fenixgroup.it

Branch Office Romania

Endosfera s.r.l. Romania - Bucarest - romanagement@fenixgroup.it

For other countries look for distributors on www.endospheres.com

ISO 13485-2011: SDS.AR.SMK.00763-14 - ISO 9001-2011: SDS.AR.SMK.00762-14

Endospheres®, Microvibrazione Compressiva®, EVA® sind eingetragene Marken der Fenix Group s.r.l.

*The results of the EVA by Endospheres are temporary and require maintenance treatments.